



This fact sheet gives you important information about the general use and care of a bathboard. It describes how to install, use and care for a bathboard.

This fact sheet does not replace specific instructions given to you by the equipment prescriber.



A bathboard can provide a safe way for a person to get in and out of a shower over the bath, however they do not fit all baths. There may be other assistive technology (equipment) available if your bath is not suitable.

A bathboard may be fitted if the bath meets the following requirements:




1. There is a flat surface on each side that the bathboard will rest on, ideally at least 30mm. (A in diagram).
2. The internal bath width (B in diagram) allows the bathboard feet to fit securely when adjusted. Make sure that the height-adjusting buttons on the legs are fully 'popped out' before sitting on the shower chair/stool.

To fit bathboard

1. Loosen feet by turning anti-clockwise.
2. Place bathboard on the bath (handle next to wall, if fitted).
3. Tighten feet by turning clockwise.
4. Test that the bathboard is secure before using.




Please note:

- Weight limit 150kg
- Temperature limit between -3° and 70°c
- Clean bathboard with warm water and a mild soap



To use bathboard

1. Stand at edge of bath with your back towards bathboard. You should feel the bath with the back of your legs.
2. Sit on bathboard as far back as you can.
3. Lift your legs into bath one at a time, holding the handle (if fitted).



Repairs to your equipment

If your equipment is broken, damaged, rusted or worn contact your local TasEquip warehouse to arrange for repair or replacement. Please call 1300 827 378 or contact your regional warehouse directly:

North

53 Dowling Street
(entry off Henry Street)
Launceston

Phone: 6777 4336

Email:

NthTasEquip@ths.tas.gov.au

North West

Mersey Community Hospital
Moriarty Road
Latrobe

Phone: 6478 5599

Email:

NWTasEquip@ths.tas.gov.au

South

Repat Centre
90 Davey Street
Hobart

Phone: 6166 7393

Email:

SouthTasEquip@ths.tas.gov.au

Review date: September 2025